

Use with textbook pages 11-14.

Microscopes

coarse focus knob	magnification power
compound-light microscope	objective lenses
electron micrograph	resolving power
eyepiece	reversed
fine-focus knob	right slide up
light source	upside-down

Use the terms in the vocabulary box to fill in the blanks. Use each term only once. You will not need to use all the terms.

- The Compound light is the microscope usually used in science classes and medical laboratories.
- The eye piece is used for viewing and contains a lens that magnifies.
- The coarse adjustment brings an object into focus at low or medium power.
- The fine adjustment brings an object into focus at high power.
- The objective have different magnification power to magnify the object.
- The light source supplies the light needed to view the slide.
- The ability to distinguish between objects that are very close together is called resolving power.
- When you look through a microscope, you will observe an image that is magnified, reversed and upside-down.
- An electron micrograph is a picture taken by a camera hooked up to an electron microscope.

Use with textbook pages 8-21.

Observing living things

Match each Term on the left with the best Descriptor on the right. Each descriptor may be used only once.

1. <u>D</u> compound microscope	A. power of the objective lens multiplied by the power of the eyepiece
2. <u>F</u> objective lenses	B. has two sets of lenses
3. <u>C</u> eyepiece	C. used for viewing and magnifying the image
4. <u>S</u> stage	D. ability to distinguish between two objects that are very close to each other
5. <u>E</u> arm	E. supports the eyepiece
6. <u>A</u> total magnification	F. have different magnification powers to magnify the object
7. <u>D</u> resolving power	G. supports the slide

Circle the letter of the best answer.

- Which of the following is not a characteristic of living things?
 - A. needs energy
 - B. hunts for food
 - C. grows
 - D. reproduces
- What is the difference between a unicellular and a multicellular organism?
 - A. size of cells
 - B. structure of cells
 - C. shape of cells
 - D. number of cells

10. Which of the following is not a use of a microscope?

- A. magnifies distant objects
- B. magnifies objects that are close together
- C. magnifies unicellular organisms
- D. magnifies bells

11. Which of the following best describes an electron micrograph?

- A. a type of electron microscope
- B. a camera hooked up to an electron microscope
- C. a picture shown on a screen hooked up to an electron microscope
- D. a scanning electron microscope

12. If the objective lens is 40x and the eyepiece lens is 10x, what is the total magnification?

- A. 400x
- B. 40x
- C. 10x
- D. 4x

13. If the objective lens is 10x and the eyepiece lens is 10x, what is the total magnification?

- A. 10x
- B. 100x
- C. 1000x
- D. 20x

Name _____ Date _____

CHAPTER 1 Parts of a Compound Light Microscope
BLM 1-5


Goal • Use this page to review the function of each part of the compound light microscope.

What to Do
Each part of the compound light microscope is listed in the left column of the table below. In the right column, describe the function of each microscope part. For assistance, refer to pages 12–13 of *BC Science 8*.

Microscope part	Function
Eye piece	used for viewing, contains a lens that magnifies
Arm	supports eye piece
Revolving nosepiece	holds 3 objectives
Objective lenses	magnify image
Fine focus knob	brings object into focus under high power
Coarse focus knob	low-med.
Stage	supports slide
Base	supports entire microscope.
Light source	supplies light.

The compound light microscope

Identify the following parts of a compound light microscope.



Key

DATE:

NAME:

CLASS:

CHAPTER 1

Estimating the Size of Microscopic Objects

BLM 1-9

Goal • Use this page to help you develop your skill of estimating the size of objects under the microscope.

Think About It

How do you estimate the size of the object you are viewing?

What to Do

Read the information below and answer the questions.

Part A: Estimating Object Size

1. Look at the four circles below. Assume that each circle below has a diameter of 2.5 cm. (Diameter is the distance across a circle.) Use the following formula to calculate the exact size of one happy face in each of the circles:

$$\text{Size of one happy face} = \text{Diameter of circle} \div \text{Number of happy faces}$$



2.5 cm



1.25 cm




0.83 cm



0.625 cm

Part B: Estimating Size Under the Microscope

- Once you know the diameter of the field of view of a microscope, you can estimate the size of the object you are viewing. The **field of view** is what you see when you look through the microscope. To find the **diameter of the field of view**, use a ruler to measure the distance across its centre. The diagram below represents a field of view when looking at millimetre markings on a ruler.



The diameter of the field of view represented on the left is 2.5 mm.

- Most objects under the microscope are much smaller than a millimetre. Try using a smaller unit, the **micrometre** (μm). Multiply the field diameter by 1000 to convert it from millimetres (mm) to **micrometres** (μm).

Convert the field of view represented above (2.5 mm) to micrometres:

The diameter of the field of view is 2500 μm .

$$1 \text{ mm} = 1000 \mu\text{m}$$

Key

DATE:

NAME:

CLASS:

CHAPTER 1

Cell Size

BLM 1-10

Goal • Use this page to practise estimating the size of cells in a field of view.

What to Do

Read the information given for each question. Answer the questions in the space provided.

1. As scientists, we must determine how small cells really are. To do this, we need to measure the diameter of the field of view.

(a) What is a field of view?

The area you see when looking through microscope

(b) What is a diameter?

The line directly through center of circle


2. When Molly looks under a microscope, before placing her specimen on the stage, she observes an empty field of view.

(a) Use your ruler to draw in the diameter of the field of view; that is, draw a line that cuts the circle exactly in half.

(b) What is the measurement of the circle's diameter?

in centimetres 1.7cm

in millimetres 17mm



3. Imagine that 10 cells of equal size fit across the diameter of the circle below.

(a) Measure the diameter of the circle. 17mm

(b) What is the span of the 10 cells? 17mm

(c) What is the span of one cell? 1.7mm

(d) Explain how you arrived at your answer for question (c).

Size of cell = Diameter / # Cells = 17mm / 10 cells = 1.7mm



4. If 10 equally sized cells fit across a field diameter of 40 mm, what is the length of one cell? Show your work.

Size = Diameter / # Cells across = 40mm / 10 cells = 4mm

